

## Available Kits for Fragment Analyzer™5200

PART NUMBER	DESCRIPTION	CONCENTRATION RANGE	SIZING RANGE
<b>RNA-Kits</b>			
<i>DNF-470-0275</i>	<i>Small RNA Kit</i>	<i>25 - 2,500 pg/μl</i>	<i>15 - 200 nt</i>
<i>DNF-472-0500</i>	<i>HS RNA-Kit (15NT)</i>	<i>50 - 5,000 pg/μl</i>	<i>200 - 6,000 nt</i>
<i>DNF-471-0500</i>	<i>RNA Kit (15NT)</i>	<i>5 - 500 ng/μl</i>	<i>200 - 6,000 nt</i>
<b>DNA-Kits</b>			
<i>DNF-474-0500</i>	<i>HS NGS Fragment Kit</i>	<i>50 - 5,000 pg/μl (smear) 5 - 500 pg/μl (fragment)</i>	<i>1 - 6,000 bp</i>
<i>DNF-473-0500</i>	<i>NGS Fragment Kit</i>	<i>5 - 100 ng/μL (smear) 0.1 - 10 ng/μl (fragment)</i>	<i>1 - 6,000 bp</i>
<i>DNF-476-0500</i>	<i>Small Fragment Kit</i>	<i>50 - 5,000 pg/μl (smear) 5 - 500 pg/μl (fragment)</i>	<i>50 - 1,500 bp</i>
<i>DNF-477-0500</i>	<i>HS Small Fragment Kit</i>	<i>5 - 100 ng/μL (smear) 0.1 - 10 ng/μl (fragment)</i>	<i>50 - 1,500 bp</i>
<i>DNF-910-K1000CP</i>	<i>CRISPR Kit</i>		<i>100 - 6,000 bp</i>
<i>DNF-467-0500</i>	<i>Genomic DNA</i>	<i>25 - 250 ng/μl (fragment)</i>	<i>75 - 60,000 bp</i>
<i>DNF-468-0500</i>	<i>HS Genomic DNA Kit</i>	<i>0.3 - 12 ng/μl (fragment)</i>	<i>75 - 60,000 bp</i>
<i>DNF-492-0500</i>	<i>Large Fragment Kit</i>	<i>0.1 - 10 ng/μl (smear) 5 - 100 ng/μl (fragment)</i>	<i>50 - 20,000 bp</i>
<i>DNF-464-0500</i>	<i>HS Large Fragment Kit</i>	<i>5 - 600 pg/μl (smear) 50 - 5,000 pg/μl (fragment)</i>	<i>75 - 48,500 bp</i>